



Reno Seat Riser

Model P325/RS | Power Chair

Reno Seat Riser

Model P325/RS | Power Chair



Controller can be fitted for left or right hand.



York Road
Bridgend Industrial Estate
Bridgend
CF31 3TB

Tel: 01656 674488
Fax: 01656 674499
Email: sales@romamedical.co.uk

www.romamedical.co.uk



We are social. Follow us on Facebook
/RomaMedicalAids

*Lower limb amputees should refer to the user manual before use.

** Passenger weight, ground conditions and battery condition/size can affect distance of travel.

*** Measurement of underside of chassis.

Roma Medical reserve the right to change specifications without prior notice.

Special Features

- 1. Luxury padded Captain Seat.
- 2. Easily dismantles without tools.
- 3. Compact and lightweight design.
- 4. Dynamic Shark controller.
- 5. On and off-board charger.
- 6. Excellent manoeuvrability.
- 7. Intelligent, regenerative, electromagnetic brakes.
- 8. Large front castor and rear wheels.
- 9. Non-marking puncture proof tyres.
- 10. Removable battery pack.
- 11. Rear mud guards.
- 12. Adjustable Height (Floor To Seat) 540-820mm / 21-32".

Specifications

Model No.*:	P325/RS
Class (Invalid Carriage):	2
Overall Length:	88cm
Overall Width:	60cm
Overall Height:	113-141cm
Seat Width:	46cm
Safe Working Slope:	6°
Maximum User Weight:	136kg
Battery Size:	22Ah
Range**:	12.5 Miles
Overall Weight:	79.5kg
Maximum Speed:	4 mph
Wheel Size - Front/Rear:	15cm / 25cm
Turning Radius:	55cm
Battery Charger:	2 Amp off board
Ground Clearance***:	5cm
Colours:	Flame Blue, Flame Red

Dealer Stamp

